


# Mitochondrial Genome Sequencing for variants identification in cancer




Group 2 - Pink

Elisabetta Callegaro  
Alice Casata  
Alessio Gambella  
Annalisa Xamin

# OUR PROPOSAL


- Efficient organelle isolation
- mtDNA extraction and recovery
- Quality Control (size, amount and purity)
- PacBio Sequel IIe NNGS (mtDNA denaturation + custom-made primer set)

No library? – No problem!



- Fractionated Mitochondria Magnetic Separation (FMMS)

- Lysis Buffer
- Phenol-Chloroform Extraction
- RNases and Proteases Treatment




- Quality Control (low-voltage agarose gel electrophoresis + Qubit + Nanodrop)



- Sequencing on PacBio Sequel IIe platform - Circular Consensus Mode (CCS)





DNA extraction



mtDNA amplification

Quality control



DNA library  
preparation


Library purification  
and  
quality control



MiSeq Illumina  
sequencing


# WetLab Workflow

1. **DNA extraction** (both nuclear and mitochondrial)
2. **mtDNA amplification**
  - a. MITO amplicons
    - i. gel electrophoresis
  - b. ROI amplicons
    - i. purification via Agencourt AMPure XP
    - ii. capillary electrophoresis via LabChip GX
3. **Library preparation**
  - a. via Nextera XT kit (MITO)
  - b. addition of adapters via PCR (ROI)
4. **Library purification and quality control**
  - a. purification via Agencourt AMPure XP
  - b. Qubit
  - c. capillary electrophoresis via LabChip GX
5. **NGS via Illumina MiSeq**




Gel electrophoresis of MITO amplicons. The MITO3 band is not visible.

1. **DNA extraction** (both nuclear and mitochondrial)
2. **mtDNA amplification**
  - a. MITO amplicons
    - i. gel electrophoresis
  - b. ROI amplicons
    - i. purification via Agencourt AMPure XP
    - ii. capillary electrophoresis via LabChip GX
3. **Library preparation**
  - a. via Nextera XT kit (MITO)
  - b. addition of adapters via PCR (ROI)
4. **Library purification and quality control**
  - a. purification via Agencourt AMPure XP
  - b. Qubit
  - c. capillary electrophoresis via LabChip GX
5. **NGS via Illumina MiSeq**



Results of capillary electrophoresis of the ROI amplicons. The control bands (without DNA template) are shorter in length and not purified.

1. **DNA extraction** (both nuclear and mitochondrial)
2. **mtDNA amplification**
  - a. MITO amplicons
    - i. gel electrophoresis
  - b. ROI amplicons
    - i. purification via Agencourt AMPure XP
    - ii. capillary electrophoresis via LabChip GX
3. **Library preparation**
  - a. via Nextera XT kit (MITO)
  - b. addition of adapters via PCR (ROI)
4. **Library purification and quality control**
  - a. purification via Agencourt AMPure XP
  - b. Qubit
  - c. capillary electrophoresis via LabChip GX
5. **NGS via Illumina MiSeq**



Results of capillary electrophoresis on Revvity LabChip GX of the MITO library.

1. **DNA extraction** (both nuclear and mitochondrial)
2. **mtDNA amplification**
  - a. MITO amplicons
    - i. gel electrophoresis
  - b. ROI amplicons
    - i. purification via Agencourt AMPure XP
    - ii. capillary electrophoresis via LabChip GX
3. **Library preparation**
  - a. via Nextera XT kit (MITO)
  - b. addition of adapters via PCR (ROI)
4. **Library purification and quality control**
  - a. purification via Agencourt AMPure XP
  - b. Qubit
  - c. capillary electrophoresis via LabChip GX
5. **NGS via Illumina MiSeq**

# Computational Analysis Workflow

1. **Reads Quality Control**
  - a. before and after adapter removal (*FastQC*)
  - b. Adapters removal (*Trimmomatic*)
2. **Alignment** of the reads with a reference genome (*BWA*)
3. **Variant calling** (*Mutect2*)
  - a. Filtering of variants (*FilterMutectCalls*)
4. **Variant functional annotation**
5. **Comparison of samples** (*PCA*)


# Results

## VEP RPE-ROI

Variant Effector Predictor - *Ensembl*

- Processed variants: 6
- Novel variants: 1
- Existing variants: 5

100% missense variants (coding)





# Results – VEP RPE–ROI

- VEP (Variant Effector Predictor) – *Ensembl*
- Considered predictors:
  - **SIFT**: based on sequence homology and the physical properties of amino acids
  - **Polyphen**: based on phylogenetic considerations

Location	Consequence	Symbol and Gene	Existing variant	APPRIS	SIFT	PolyPhen	Clinical Significance
M:3739-3739	Missense variant	MT-ND1 ENSG00000198888	-	P1	0.23	0.119	-

# Results


## VEP RPE-MITO

### Variant Effector Predictor – *Ensembl*

- Processed variants: 21
- Novel variants: 1
- Existing variants: 20

50% missense variants


50% synonymous variants



# Results – VEP RPE-MITO

Location	Consequence	Symbol and Gene	Existing variant	APPRIS	SIFT	PolyPhen	Clinical Significance
M:13145-13145	Missense variant	MT-ND5 ENSG00000198786	rs386829175	P1	1	0.003	benign
M:6419-6419	Missense variant	MT-CO1 ENSG00000198804	rs1603220461	P1	0	0.999	-
M:9912-9912	Missense variant	MT-CO3 ENSG00000198938	rs28580363	P1	0	0.957	-
M:15326-15326	Missense variant	MT-CYB ENSG00000198727	rs2853508	P1	0.21	0.009	benign, likely pathogenic
M:8860-8860	Missense variant	MT-ATP6 ENSG00000198899	rs2001031	P1	0.43	0.003	benign

# Results – PCA



## Expectations

- ★ Healthy and cancer cell lines are separated
- ★ MITO and ROI cluster together

## Limitations

- ★ Overlapping points
- ★ Random distribution (Jitter function)

## Explanations

- ★ Similar data
- ★ Problems in the PCA

# VCF file analysis

**Table S1** | Mutational burden of cancer vs healthy cell lines. The variations have been counted.

	A549	Cal51	HCT116	Nalm6	RPE
MITO	36	26	34	57	27
ROI	8	10	7	21	8
MITO $\cap$ ROI	39	32	38	67	30

## Conclusions

- ★ Mutational burden higher in cancer cell lines
- ★ ROIs were valuable to identify variants

# VCF file analysis

Only one variant is present in all tumor cell lines but not in RPE cells



**MT-NDR2 gene**

chrM:2.465

Total count: 5024

A : 760 (15%, 33+, 727- )

C : 4 (0%, 4+, 0- )

G : 19 (0%, 7+, 12- )

T : 4241 (84%, 2949+, 1292- )

N : 0

-----

DEL: 8

INS: 1

Table S2 | Position and type of the SNPs found in RPE samples but not in tumor cells.

Position in ChrM	Gene	Type	RPE MITO	RPE ROI
152	-	Upstream/downstream gene variant	yes	yes
1959	MT-RNR2	Non-coding transcript exon variant	yes	no
3363	MT-ND1	-	no	yes
5318	MT-ND2	Synonymous variant	yes	no
5691	MT-TN	-	no	yes
9912	MT-CO3	Missense variant	yes	no
9950	MT-CO3	Synonymous variant	yes	no
13145	MT-ND5	Missense variant	yes	no
15466	MT-CYB	Synonymous variant	yes	no
15721	MT-CYB	Synonymous variant	yes	no
16192	-	Upstream/downstream gene variant	yes	no

# OUR PROPOSAL

## VS

# THE HARSH

# REALITY

## *Main Divergences*



### Distinct sequencing technologies:

- *Illumina MiSeq*
- *PacBio Sequel IIe*

### Distinct extraction processes:

- *Nuclear/mitochondrial DNA extraction + enrichment*
- *Organelle isolation + mtDNA extraction*

### Distinct quality control approaches and timing

### Library preparation



## *Some Commonalities*

### Quality control techniques:

- *Qubit fluorometer*
- *Gel electrophoresis*

# CRITICALITIES AND PROBLEMS

*(THE DAILY BREAD OF ANY  
RESPECTABLE SCIENTIST)*

## *Notable Mentions*

- *Unsuccessful MITO3 amplification*
- *Unsatisfactory library size and purity*
- *Irrelevant PCA plot*





*THANK YOU  
FOR  
YOUR  
ATTENTION!*



# TEAM



# PINK